



CLA 9/7

Scopin-Broncho-Boy II

by Prof. Dr. med. John A. Nakhosteen.

CLA 9/7

Scopin-Broncho-Boy II

Universal Bronchoscopy and Anesthesiology Trainer produced in co-operation with Prof. Dr. J. A. Nakhosteen. Life-sized, in plastic.

While retaining the life-like anatomical features of the original, the Broncho-Boy II incorporates a number of innovations enabling its use in a variety of training settings, including:

Components:


1. SCOPIN Broncho Boy II CLA 9/7 with head in movable assembly.
2. Warning signal emitter in upper jaw, nasopharynx, tracheobronchial tree with twist lock.
3. Trunk with elastic supports and auscultation membranes, mounted on wooden base with folding hinge.
4. Removable healthy and fluorescent bronchial tree, with twist lock.
5. Removable pulmonary alveoli with membrane connection with twist lock.
6. Alarm tone generator.
7. Intubation tube
8. Silicone oil, 250 ml pressure spray bottle
9. Instructions for use



CLA 9/7, opened

Each CLA SCOPIN model receives a serial number, which may be found inside the thorax on the right side. Please always quote this number in full when re-ordering or in case of complaint.

CLA 9/7
608.09.19/FR



Serial number in the chest (right inside)

CLA reserves the right to modify models as required by technical and scientific developments.

Special equipment on request:

Aluminium transport and storage case CLA 9/2-Z



Construction

A specially developed casting process has enabled exact reproduction of internal structures and contours of mucous membranes. Anatomical landmarks, this cast, aid the trainee bronchoscopist in mastering problems of endoscopic orientation. For those experienced with the rigid bronchoscope and wishing to learn flexible endoscopy, practicing with the SCOPIN model will greatly ease the switch.


The following characteristics typify this unique training model:

1. Retroflexion and left-right rotation of the head for laryngoscopic, oro-tracheal and bronchoscopic intubation.
2. Acoustic warning signal if excessive pressure exerted on upper teeth by rigid tube or laryngoscope
3. Slightly widened and softer nasal passages for fibre-optic guided, naso-tracheal intubation; reduced danger of damage to sheath.
4. Interchangeable, membrane-connected bladders (Broncho-Boy II AN) enabling stethoscopic confirmation of ET-tube placement.
5. New, spring-based supports for tracheo-bronchial tree, increasing pliability, reducing likelihood of bronchoscope damage.
6. Exchangeable fluorescent tracheo-bronchial system with »cold« areas of reduced fluorescence for practising autofluorescent bronchoscopy.

Adjusting incisor acoustic trigger device

Excessive force on the upper incisors by laryngoscope or rigid bronchoscope activates an **acoustic signal**, the trigger-sensitivity of which is adjustable. The adjusting screw is located next to the right upper molars. **Turning to the right**, i.e., tightening, **reduces** sensitivity. **Turning to the left**, i.e., loosening the screw, **increases sensitivity**. (A second screw in the center of the palate secures the upper jaw).

The housing for the signal-device, located within the thorax at the lower right, also contains an **alkaline battery** with a service life of about one year. To replace the battery, push the entire housing downward, remove it and disconnect the two snap connections. This gives access to the battery, which can now be replaced.



Adjusting the sensitivity of the switch, which triggers off the acoustic signal.



Alarm tone generator, dismantled

Intubating the Broncho-fiberscope

It is recommended that beginners spend 10 hours practising with the Nakhosteen Bronchoscopy Model.

Preparing for use

1. Removal of sternal window: **begin** at lower end
2. Undo the transportation lock and turn knurled nut upwards, as far as it will go.
3. Removing the Silicon lubricant aerosol
4. Spraying the oral and pharyngeal area
5. Spraying the nasal passage
6. Connecting the healthy or the fluorescent lower part of the bronchial tree with the upper part
7. Replacing sternal window: begin with **upper** end

Before each intubation, the interior of the nose, pharynx and the fiberscope are to be slightly moistened with silicone oil. Never use paraffin spray!

A semi-reclining position is preferred so that the SCOPIN is facing the examiner. The SCOPIN is laid flat in order to enable the examiner to intubate from the head end of the model. If the student is already familiar with the rigid bronchoscope, placement from the head end of the model will pose less difficulty.

There are basically two types of fiberscope intubation:

trans-nasal and trans-oral. Trans-oral access should be performed by way of a bite block or intubation tube. The student should become acquainted with both methods. Trans-nasal intubation is more comfortable for the patient, so if there is no specific indication for it to be conducted by the oral route, nasal insertion is preferred.

The intubation may proceed either transnasally or transorally, and the transoral approach must be either through a gag or ring (to prevent damage through biting), or through a flexible endotracheal tube. The trainee should gain proficiency in all these approaches. As far as the patient is concerned, however, the transnasal route is best tolerated, so unless there is a special indication for an alternative transnasal intubation should be used.

The smallest diameter of both nasal passages of the Nakhosteen Bronchoscopy Model is 6,5 mm, so that every currently available bronchofiberscope can be passed through the nose. In such narrow spaces, however, the widest passage must be sought and found by careful manipulation of the scope. In principle, passage of the instrument should be possible either above or below the inferior turbinate bones, but it may well be - as is often the case with patients - that one of these routes is better. These variations formed a major part of the deve-

lopmental concept of the SCOPIN model.

Since light can enetrate the SCOPIN tracheobronchial tree, the position of the endoscope can be checked at will, simply be removing the sternum. Please note that in opening the sternal window the lower end is removed first, and in closing the window the lower end is replaced last.

Access to the lingula bronchus was purposely made somewhat more difficult. Firstly, in patients the lingula is sometimes quite difficult to reach; and secondly, practicing here may be considered useful for the mire advanced trainee.

Some general hints for training in bronchofiberscopy using the Nakhosteen Bronchoscopy Model, SCOPIN

Guidance

It is important that the trainee practicing with the SCOPIN Model have a thorough knowledge of the more important anatomical landmarks of the nose, throat, and the tracheo-bronchial tree.


Seen from the ventral aspect, the epiglottis is "like a small tongue stretched toward the examiner".

The vocal cords are like "the tip of an arrow pointing ventrally".


The fibrous membrane (pars membranacea) coursing longitudinally in the tracheo-bronchial system is **dorsal**, and the horse-shoe shaped cartilagenous tracheal rings are **ventral** and **lateral**.

Both the hand and the eye must be conditioned in this training process: the scope, and the eye in recognizing anatomical landmarks. The two hand movements constantly interacting with one another are axial rotation and bending and extending the tip of the instrument.


Use extreme care in maneuvering the flexible bronchoscope; any forceful manipulation may lead to damage and costly repairs.




ill. 1
Vocal cords and ventral commissure



ill. 2
Trachea; main bifurcation distally



ill. 3
Right upper lobe bronchus; subsegments of RB 3 visible



ill. 4
Main stem bronchi; branchings to right upper lobe and middle lobe bronchi visible


Care

The inner and outer surface of the SCOPIN-Model requires no special care. When dirty, simply clean with water or soap solution.

Treatment of the nasal, oral and pharyngeal cavity, as well as the oesophagus and bronchial tree, must be carried out in the same way. However, after each cleaning, it is desirable to spray on a very fine layer of silicone oil, which is included, to achieve better insertion.


Tape residues are removed with petroleum ether.

If markings are to be made, use only pencil. Traces of crayons, pens, markers, and other coloured solutions, cannot be removed.



**before
intubation
always...**

**Fully unscrew the
two nuts inside the
trainer, so what...**



**... The head can be
retroflexed. And ...**



**... Coat instrument
surfaces with sili-
cone oil or liquid
prior to intubation**

Laryngoscopic Intubation



Laryngoscope is used to visualize the vocal cords, then the bronchoscope intubated


Direct Intubation



The epiglottis, visualized in the mid-line, is gently levered upward with the tip of bronchoscope, examiner's thumb acting as fulcrum




Once vocal cords come into view, the bronchoscope is carefully advanced into trachea




Intubated open-tube bronchoscope, trainer's head can be rotated 45° left or right as needed

Flexible Fiberoptic Bronchoscopy (FFB)




Trans-nasal Flexible Fiberoptic Bronchoscopy



Trans-oral Flexible Fiberoptic Bronchoscopy




Fluorescent tracheal bronchial tree and exchangeable lung bags with membrane connection



Exchanging the auscultation bellows with tracheo-bronchial tree "Dot do Dot" 90° clock-wise turn

Laryngoscopic Intubation



Laryngoscopic intubation of endotracheal tube



Auscultation following laryngoscopic intubation

Flexible Fiberoptic-guided Intubation




Flexible fiberoptic-guided naso-tracheal intubation



Flexible fiberoptic-guided oro-tracheal intubation

RIGHT

LEFT



Accessories and spare parts list for CLA 9/7

CLA 97100

Head with upper part of the bronchial tree, exchangeable eyes, set of teeth with signaling device, cable with connecting plug for tone generator, 2 fastening screws with self-locking nut

CLA 97110

1 pair of eyes with lids

CLA 97111

Eyelids (1 pair)

CLA 97112

Eyelids, on the right

CLA 97113

Eyelids, on the left

CLA 97114

Glass eyes (1 pair)

CLA 97120

Healthy bronchial tree, lower part

CLA 97140

Fastening screws (1 pair) with lock nuts M 6

CLA 97200

Thorax with flexible neck, chest wall, fixing hinge, support for the bronchial tree, tone generator plate, 2 membranes for auscultation areas, cable clips and screwdriver retainer

CLA 97210

Support for the bronchial tree

CLA 97230

Chest cover

CLA 97240

Sound generator with housing, connection jack and battery

CLA 97262

Screwdriver 2,5 mm

CLA 97270

Base plate (670 x 380 x 10 mm) with 2 guide bars, knurled nut and drill holes for the fixing hinge, with rounded corners and bevelled edges

CLA 97300

Bronchial tree connector piece with lungs and connecting tubes for intubation and ventilation

CLA 97310

Bronchial tree connector piece

CLA 97320

Lung bags, 1 pair

CLA 97321

Lung bag, right

CLA 97322

Lung bag, left

CLA 97330

2 connection tubes from the bronchial branch to the membranes

CLA 97700

Operations instructions, german

CLA 97702

Operations instructions, english

CLA 97704

Operations instructions, french

CLA 97708

Operations instructions, spanish

CLA 800/1

Silicon separating spray, 250 ml net contents (pressure spray bottle)

CLA 9/8

Fluorescing Tracheo-Bronchial Tree with regions of reduced Fluorescence

Special equipment on request:

CLA 9/2-Z

Aluminium transport case with foam rubber
Height: 28 cm, width: 76 cm, depth: 44 cm, weight: 6,0 kg.

CLA 97130

Sick bronchial tree, lower part



CLA 9 /7
SCOPIN-Broncho-Boy II
 disassembled, with a note containing order numbers
 for accessories and spare parts

**Special equipment
on request:**



CLA 9/2-Z

Aluminium transport case
with foam rubber
Height: 28 cm, width:
76 cm, depth: 44 cm,
weight: 6,0 kg.




CLA 97130

Sick bronchial tree, lower
part



Detail from CLA 97130

Tumor occluding bronchus intermedius; right upper lobe displays substance for biopsy (not shown).



Detail from CLA 97130

Adenoma occluding left lower lobe bronchus. Left upper lobe bronchus partially, LB 6 wholly, visible.

CLA 9/8

Fluorescing Tracheo-Bronchial Tree with regions of reduced Fluorescence

Produced in co-operation with Prof. Dr. J. A. Nakhosteen,

In keeping with the 20-year-old tradition of up-dating the Nakhosteen Broncho Boy Bronchoscopy Teaching Model, the latest innovation is the auto-fluorescing endoscopic systems for use with the Light Imaging Fluorescence Endoscope (LIFE® Xillix/Olympus), or SAFE (1000 Pentax).



CLA 9/8

The bronchial mucosa appears normal with conventional (white light) bronchoscopy, but on being illuminated in the fluorescent mode, emits a green image identical to that of normal mucosa.

The trachea and left-sided bronchial branches appear normal in the fluorescence mode, but distinct areas of reduced fluorescence can be observed on the right-sided bronchial tree (upper lobe spur, middle lobe spur, right main bronchus and RB 9/10 spur).


Endoscopists using the LIFE-System or SAFE 1000 for the first time can practice switching from white light to fluorescent mode, and in the latter mode, learn <gain> (light intensity) adjustment, in addition to recognizing areas of reduced fluorescence.

The Fluorescent Tracheo-Bronchial Tree can be added to each CLA Broncho Boy with the order numbers CLA 9 to CLA 9/7 or can be supplied additionally as the interchangeable lower part of the tracheo-bronchial tree. Easy to interchange, thanks to bayonet fixing



III. 5

White light appearance of fluorescing bronchial tree: RB 9/10 spur with normal features



III. 6

Fluorescence bronchoscopic appearance of same spur as in Fig. 5 (RB 9/10): in the middle is a cold spot corresponding to a Carcinoma in situ, visible only in fluorescence mode